

# FT-IR SAMPLING TECHNIQUES

Shorter Wavelengths

Near-IR

Mid-IR

Far-IR

Longer Wavelengths

UV-Vis

## X-H Overtones and Combination Bands

## Fundamental Vibrations, Overtone and Combination Bands


## Fundamental and Heavy Atom Vibration Bands

Microwave

Wavelength  
0.7 μm  
Wavenumber  
14000 cm<sup>-1</sup>  
Frequency  
10<sup>14</sup> MHz  
Energy  
1.771 eV

### Transmission/Absorption

Sample is positioned directly in the beam of the FT-IR spectrometer



Beer-Lambert Law  
 $A = a \cdot b \cdot c$

A - Absorbance  
a - absorptivity  
b - pathlength  
c - concentration

Absorbance is linearly proportional to sample concentration and pathlength.

#### Attributes


- High quality, representative spectral data, compatible with many digital libraries
- Spectra represent bulk sample composition
- Ability to control sample thickness – intensity of spectral features/sensitivity

#### Special Considerations

- Thick samples may absorb totally and need to be thinned-out or diluted
- Original material may be destroyed in sample preparation process
- Pellet making and other sample preparation methods require skill and time
- Beer's law is limited to simple, linearly absorbing samples

#### Applications

- Solids, gels, pastes, liquids and gases
- Film thickness and liquid cell pathlength
- Qualitative and quantitative



Transmission Ranges of Popular IR Windows

NaCl	52600-455 cm <sup>-1</sup>
KBr	48000-345 cm <sup>-1</sup>
CaF <sub>2</sub>	79500-895 cm <sup>-1</sup>
BaF <sub>2</sub>	66600-690 cm <sup>-1</sup>
ZnSe	15000-460 cm <sup>-1</sup>

Free Standing Film Thickness (μ)  
 $T = (N \cdot 10,000) / (2n \cdot \Delta\nu)$

T - thickness of polymer film (μ)  
Δν - wavenumber range of fringes count (cm<sup>-1</sup>)  
N - number of fringes within Δν  
n - refractive index of polymer

Cell Pathlength (mm)  
 $P = (N \cdot 10) / (2 \cdot \Delta\nu)$

P - cell pathlength (mm)  
Δν - wavenumber range of fringes count (cm<sup>-1</sup>)  
N - number of fringes within Δν  
n - refractive index of polymer


#### Sample Handling

- Films, thin materials – direct placement in the transmission sample holders
- Powder and solid samples: KBr pellets, (solid samples require processing), mulls
- Gels/pastes: smears on IR transparent windows
- Liquids: liquid cells
- Gases: short- and long-path gas cells

2.5 μm  
4000 cm<sup>-1</sup>  
10<sup>14</sup> MHz  
0.496 eV

### Attenuated Total Reflection (ATR) (internal reflection)

Sample is placed in intimate contact with the internal reflection element (IRE)




Depth of Penetration  
 $d_p = \frac{\lambda}{2\pi(n_1^2 \sin^2 \theta - n_2^2)^{1/2}}$

d<sub>p</sub> - depth of penetration  
λ - wavelength of light  
θ - angle of incidence  
n<sub>1</sub> - refractive index of crystal  
n<sub>2</sub> - refractive index of sample

Depth of penetration can be controlled by changing angle of incidence beam and/or use of crystals with different refractive indices.

#### Number of Reflections



$$r_s = L / (2T \tan \theta)$$

r<sub>s</sub> - number of reflections  
L - crystal length  
T - crystal thickness  
θ - angle of incidence

#### Refractive Index and Critical Angle



$$\theta_c = \sin^{-1} \left( \frac{n_2}{n_1} \right)$$


θ<sub>c</sub> - critical angle  
n<sub>1</sub> - refractive index of crystal  
n<sub>2</sub> - refractive index of sample

#### Attributes

- Minimal sample preparation
- Short sampling time and easy cleaning
- High versatility, non destructive sampling

#### Special Considerations

- The angle of incidence must exceed the critical angle to produce good quality spectra
- ATR is not a bulk technique; it measures only the surface of a potentially inhomogeneous sample
- Correction needs to be applied to spectra for direct comparison with transmission data
- ATR crystals must be cleaned between samples to avoid contamination and carry-over
- Spectral ranges and properties of ATR crystals vary and need to be considered when working with specific samples



Transmission Ranges of Popular ATR Crystals

ZnSe	52600-650 cm <sup>-1</sup>
Ge	5500-780 cm <sup>-1</sup>
Diamond	30000-50 cm <sup>-1</sup>
Si	8900-1460 cm <sup>-1</sup>

#### Applications

- Solids, gels, pastes and liquids
- Qualitative and quantitative


#### Sample Handling

- Intimate contact between the sample and the ATR crystal is required

25 μm  
400 cm<sup>-1</sup>  
10<sup>11</sup> MHz  
0.050 eV

### Diffuse Reflection (scattered reflection)

Sample is diluted with KBr (or similar) and placed in a cup for analysis  
Large surface areas can be analyzed directly



#### Attributes

- Simplified sample preparation
- Direct measurement of large samples
- Representative sample spectra (not a surface measurement)
- Sensitivity and versatility


Kubelka-Munk Equation

$$f(R) = \frac{(1-R)^2}{2R} = \frac{k}{s}$$

R - absolute reflectance of the layer  
k - molar absorption coefficient  
s - scattering coefficient

#### Special Considerations

- Diffuse reflection spectra will appear different to their transmission equivalents
- Kubelka-Munk conversion brings spectra to transmission-like format
- If not eliminated, specular reflection component may distort spectra
- Spectra of samples with high refractive indices require higher dilution to reduce distortion



#### Applications

- Measurements of powdered samples
- Selected liquids/liquid residues
- Rough surface solids/coated surfaces
- Abrasion measurements
- Temperature and catalytic studies



#### Sample Handling

- Powders – 1% to 3% dilution in KBr, KCl or similar matrix
- Small particle size (50 μ or less) and homogeneous mixture
- Maintain more than 2 mm sample depth, even and loose packing in the cup
- Liquids – dispense small amount on top of KBr powder located in the sampling cup
- Large samples placed directly on top of upward looking diffuse reflection accessory or integrating sphere

### Specular Reflection (mirror-like reflection)

Sample is placed on top of the specular reflection accessory

#### Reflection

- Incident beam angle = reflected beam angle
- Angle of incident beam transmitted into a sample is reflected by Snell's Law

#### Reflection-Absorption at Near-Normal and Grazing Angle (< 65°)

- Beam enters and exits sample deposited on reflective substrate
- Equivalent to double-pass transmission measurement

#### Reflection in Free Standing Films

- Beam enters sample and reflects from top and bottom surfaces of the film
- Constructive and destructive interference generate fringe pattern

I<sub>0</sub> - Incident beam  
I<sub>t</sub> - Transmitted beam  
I<sub>r</sub> - Reflected beam (1)

I<sub>2</sub> - Reflected beam (2)  
θ<sub>i</sub> - angle of incidence  
θ<sub>t</sub> - angle of trans. light

n<sub>1</sub> - refr. index of air  
n<sub>2</sub> - refr. index of sample

Snell's law indicates that the angle of the incident beam into the sample is equal to the angle of reflection:  $n_1 \sin \theta_i = n_2 \sin \theta_r$

#### Attributes

- No or limited sample preparation
- Measurement simplicity/versatility
- Sensitivity/ability to look at very thin samples
- Non-destructive sampling


Film Thickness (μm)

$$T = \frac{10,000 \cdot N}{(2 \cdot \Delta\nu)(n^2 \sin^2 \theta)^{1/2}}$$

T - thickness of film (μm)  
Δν - wavenumber range of fringes count (cm<sup>-1</sup>)  
N - number of fringes within Δν  
n - refractive index of polymer

#### Special Considerations

- Specular reflection spectra of thick reflecting samples may be distorted due to anomalous dispersion. Kramers-Kronig transform brings spectra to absorption-like format.



#### Applications

- Measurements of thin and thick film composition and thin film thickness
- Analysis of monolayers on reflective substrates
- Reflectance measurements – coatings, optical components/glass, EPI, lubricants on hard drives

#### Sample Handling

- Polarized light is recommended for grazing angle measurements

Think  
**PIKE**  
TECHNOLOGIES